

# GO-GO<sup>®</sup>

## ULTRA X

### **Economy Class Price. First Class Features.**

The **Go-Go<sup>®</sup> Ultra X** is loaded with design features like one hand Feather-touch disassembly and a convenient drop-in battery box that makes transport and travel worry free. The Go-Go<sup>®</sup> Ultra X delivers these features and more, along with the performance you expect from the first name in travel mobility, making it the ultimate travel scooter value.



*Go-Go<sup>®</sup> Ultra X, 3-wheel*



*Go-Go<sup>®</sup> Ultra X, 4-wheel*

**Available as a 3 & 4 wheel scooter**



**Feather-touch disassembly.**

**Pride<sup>®</sup>**  
Mobility Products Corp.



### Specifications

Model numbers	SC40X (3-wheel); SC44X (4-wheel)
Max. weight capacity	260 lbs.
Maximum speed <sup>1</sup>	Up to 4.25 mph
Ground clearance <sup>2</sup>	3-wheel: 2.625" 4-wheel: 2.5"
Turning radius <sup>2</sup>	3-wheel: 32.5" 4-wheel: 43.25"
Overall length <sup>2</sup>	3-wheel: 36.5" 4-wheel: 39.75"
Overall width <sup>2</sup>	3-wheel: 19.25" 4-wheel: 19.125"
<b>Tires</b>	
Front	3-W: 2" x 8" solid; 4-W: 2" x 7" solid
Rear	2.5" x 8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drive train	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 2 amp
Per-charge range <sup>1,5</sup>	Up to 9 miles (w/ 2 x 12 AH)
Battery requirements <sup>3</sup>	(2) 12V deep-cycle (standard)
Compatible batteries	12 AH
<b>Weight</b>	
Batteries (2 required)	12 AH (21 lbs. battery pack) <sup>4</sup>
Weight w/o battery pack <sup>2</sup>	3-wheel: 71 lbs. 4-wheel: 76 lbs.
Heaviest piece <sup>2</sup>	3-wheel: 28 lbs. (front section) 4-wheel: 33 lbs. (front section)

### Standard Colors



Red



Blue

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>AGM or Gel-Cell type required

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>This range was obtained with the 4-wheel model.

### Features

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black non-scuffing tires
- Auto-connecting front to rear cable
- Front frame-mounted seat post offers maximum stability
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board


## Silver Star Lifts & Ramps

These Silver Star® Lifts are compatible with the Go-Go® Ultra X.

- Backpacker® Plus
- Backpacker® MV
- Power Tote Mini (3-wheel only)



Full Platform\*



Silver Boom 300

*NOTE: All specifications subject to change without notice.*

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

**Pride**<sup>®</sup>  
Mobility Products Corp.