

Jazzy[®] SELECT

The highly maneuverable Jazzy Select[®] delivers the advantages of in-line motor technology for enhanced efficiency, torque, range and performance. Plus, it's loaded with performance and convenience features like standard Active-Trac[®] Suspension that make it as easy to enjoy as it is to use.



Jazzy Select[®] shown
in Jazzy Red with gray vinyl
high-back seat, PG GC 2
controller and foot platform



800-800-8586 (US) • 888-570-1113 (Canada)
Visit us on the web: www.pridemobility.com

Pride[®]
Mobility Products Corp.

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	2" (1.875" at center of frame)
Turning radius	20"
Overall length	35.5" without front riggings
Overall width	23.75"
Wheels	Drive: 10"
Front anti-tips	5" solid
Rear casters	6" solid
Suspension	Active-Trac®
Drive train	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 2 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	99 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

COLOR CHOICES



Jazzy Red



Jazzy Blue

* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.

800-800-8586 (US) ■ 888-570-1113 (Canada)
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

MODELS

Jazzy Select 2S-C
Jazzy Select 2S-SS

FEATURES

- 40 amp, PG GC 2 controller
- Built with ease of service in mind
- Easy side access to batteries
- Tool-less front panel access to electronics to minimize service time
- Active-Trac® Suspension

OPTIONS

- Depth-adjustable fold flat seat provides fast, simple adjustment for enhanced, personalized comfort*
- High-back seat with headrest
- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Silver Star
Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker MV
- Silver Boom 300



PRIDE
PERFORMANCE MOBILITY

— LIVE YOUR BEST™ —