

# Jazzy<sup>®</sup>

## SELECT 6 ULTRA

The Jazzy<sup>®</sup> Select<sup>™</sup> 6 Ultra provides the advantages of in-line motor technology, maximum front stability and the environmental access benefits of a standard power elevating seat. These features, along with Active-Trac<sup>®</sup> Suspension and 6" front and rear casters make the Jazzy Select 6 Ultra high performing and easy to use.



Jazzy<sup>®</sup> Select<sup>™</sup> 6 Ultra shown in Jazzy Red with contemporary design high-back seat, PG VR2 controller and foot platform


800-800-8586 (US) ■ 888-570-1113 (Canada)  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride<sup>®</sup>**  
Mobility Products Corp.

## SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	2"
Turning radius	20.5"
Overall length	36.625" without front riggings
Overall width	23.5"
Wheels	Drive: 10"
Front casters	6" solid
Rear casters	6" solid
Suspension	Active-Trac®
Drive train	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG 6-key VR2 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 15 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	24.5 lbs. each
Base weight	114 lbs.
Standard seat weight	37 lbs. (medium-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors

## COLOR CHOICES



Jazzy Red



Jazzy Blue

## MODELS

Jazzy Select 6 Ultra-C  
Jazzy Select 6 Ultra-SS

## FEATURES

- Standard power elevating seat
- 6" casters on front and rear for maximum stability
- Dual in-line motor technology for enhanced efficiency, torque and range
- Battery access through side panel without removal of the seat
- Tool-less front panel access to electronics to minimize service time

## OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## RECOMMENDED LIFTS

**Silver Star®**  
**Lifts & Ramps**

- Backpacker® MV
- Backpacker® Plus
- Full Platform
- Silver Boom 300



**Pride®**  
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada)  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.