

Stylus®



Quality and Comfort...for Life.

The Stylus® is designed to meet the needs of active users. It's extremely lightweight, highly durable, and provides enhanced fit and function. The Stylus offers all this in a very cost effective package.



PRIDE®
PERFORMANCE MOBILITY

— LIVE YOUR BEST™ —

SPECIFICATIONS

PRODUCT WEIGHT 33.5 lbs. without leg rests¹

WEIGHT CAPACITY 250 lbs.

FRAME MATERIAL Lightweight aluminum

SEAT-TO-FLOOR HEIGHT² 19" (standard)

SEAT WIDTHS 16", 18" and 20"

SEAT DEPTHS 16" and 18"

BACK HEIGHTS 17" - 19"

REAR WHEELS 24" (standard) composite mag wheels with polyurethane tires and aluminum handrims

WHEEL LOCKS Push-to-lock

CASTERS 8" flat-free

ANTI-TIP WHEELS Removable, adjustable

WARRANTY Lifetime limited on side-frames and cross-frames;
1-year limited on materials

FEATURES

- Adjustable tension seat upholstery increases comfort and positioning
- Multiple-position rear wheel and caster mounts optimize center-of-gravity and seat-to-floor height adjustments
- Quick-release, removable rear wheels
- Full-profile 24" x 1 3/8" rear wheels for enhanced outdoor capabilities
- Long handle wheel locks ensure comfort and security
- Industry-compatible leg rests with aluminum footplates for durability

OPTIONS

- Synergy[®] Cushions
- Synergy[®] Shape Back
- Angle-adjustable back canes
- Tension adjustable back upholstery
- Adjustable rear anti-tips
- Oxygen holder
- Elevating leg rests with calf pads
- Height-adjustable armrests, desk or full length
- Flip-up, desk length armrests
- 6" flat-free casters
- 22" rear wheels
- Pull-to-lock wheel locks
- Wheel lock extensions

RECOMMENDED LIFTS

Silver Star[®]
Lifts & Ramps

- Caddy
- Silver Boom 300



STYLUS
shown with Synergy[®] Solution[®] cushion

Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST™ —

(US) 800-800-8586 ▪ (Canada) 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.

¹) Standard 16" x 16" seat size and flip-up, desk-length arms

²) Height varies from 17" - 20" based on wheel size