

VICTORY[®] 9

THE ULTIMATE 9" SCOOTER
WITH LUXURY FEATURES AT
AN ECONOMY PRICE.

VICTORY[®] 9

The **Victory[®] 9** delivers the high performance operation and feature-rich design of the Victory series of scooters. Feather-touch disassembly, exclusive low-profile tires, back-lit battery gauge, and wraparound delta tiller all combine to make the Victory 9 an exceptional value.



Pride[®]
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Victory® 9	3-wheel	4-wheel
Model number:	SC609	SC709
Weight capacity:	300 lbs.	300 lbs.
Maximum speed ¹ :	Up to 5.25 mph	Up to 5.25 mph
Ground clearance ² :	2"	2.25"
Overall length ² :	43"	45.5"
Overall width ² :	22.25"	22.5"
Turning radius:	42"	51.5
Range per charge ^{1,6} :	13 miles ⁵	
Tire Size:	Front: 8" solid Rear: 9" solid	Front: 9" solid Rear: 9" solid
Total Weight Without Batteries ^{2,7} :	117.5 lbs.	131 lbs.
Heaviest piece when disassembled ² :	Front: 48.5 lbs.	Front: 59 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black vinyl Dimensions: 18" overall width 17" effective depth	
Seat-to-Deck Range:	17.25" - 18.75"	
Drive System:	Rear-wheel drive, sealed transaxle, 24-volt DC motor	
Dual Braking System:	Regenerative and electromechanical	
Battery Requirements ^{3,4} :	Two 12-volt, deep cycle Size: U-1 Weight: U-1:24.5 lbs. each	
Battery Charger:	3-amp, off-board	

Colors:



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required

⁴Battery weight may vary based on manufacturer.

⁵This range was obtained with the 4-wheel model.

⁶Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁷Includes standard seat weight

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Features

- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Pride's exclusive stylish, light-weight, low profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges & eight-position seat rotation
- 300 lbs. weight capacity

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Silver Star Lifts & Ramps ★

- Backpacker® Plus
- Full Platform
- Silver Boom 300
- Backpacker® MV

